

ROSENBAUM POCKET VISION SCREENER

95

874

2843

638 E W E X O O 14 10 $\frac{14}{20}$

8745 E M W O X O 10 7 $\frac{10}{20}$

63925 M E E X O X 8 5 $\frac{8}{20}$

428365 W E M O X O 6 3 $\frac{6}{20}$

317218 W W X X O 5 2 $\frac{5}{20}$

217218 W W X X O 4 1 $\frac{4}{20}$

117218 W W X X O 3 1+ $\frac{30}{20}$


distance
equivalent

Point
Jaeger $\frac{14}{20}$
 $\frac{10}{20}$
 $\frac{8}{20}$
 $\frac{6}{20}$
 $\frac{5}{20}$
 $\frac{4}{20}$
 $\frac{30}{20}$

Card is held in good light 14 inches from eye. Record vision for each eye separately with and without glasses. Presbyopic patients should read thru bifocal segment. Check myopes with glasses only.

DESIGN COURTESY J. D. ROSENBAUM, M.D.

PUPIL GAUGE (mm.)



Catalog No. 13-1243-1
GRAHAM-FIELD SURGICAL CO., INC.

PRINTED IN U.S.A.