

## DENTAL EQUIPMENT MAINTENENCE

DESCRIPTION OF EQUIPMENT	DAILY	WEEKLY	BI-WEEKLY	MONTHLY	BI-MONTHLY	SEMI-ANNUAL	ANNUAL	AS NEEDED
<b>Automatic Film Processor</b>								
a.)Check Chemistry								
b.)Change Chemistry								
c.)Clean Racks								
d.)Clean Entire Unit								
<b>Chairs</b>								
a.) Clean								
<b>Compressor</b>								
a.)Check Moisture Indicator								
b.)Check Air Intake Filters								
c.)Regenerate Dessicant								
d.)Check Oil Level								
e.)Change Oil								
<b>Dental Operating Light</b>								
a.)Clean Reflector								
b.)Check and Clean Sheild								
c.)Check for Drift								
<b>Nitrous System</b>								
a.)Check Pressure Overnight								
b.)Turn System On/Off								
c.)Check Rubber Goods for Wear								
<b>Sterilizer, Pelton</b>								
a.)Activate Safety Valve								
b.)Clean Chamber								
c.)Spore Test								
d.)Check Door Seals								
<b>Sterilizer, Statim</b>								
a.)Lubricate Cassette seal								
b.)Apply Stat Dry to Cassette								
c.)Replace Air Filter								
d.)Empty Used Water Container								
<b>Dental Unit</b>								
a.)Clean Out Vacuum Trap								
b.)Clean Vacuum Lines								
c.)Flush Handpiece Tubings								

## DENTAL EQUIPMENT MAINTENENCE

DESCRIPTION OF EQUIPMENT	DAILY	WEEKLY	BI-WEEKLY	MONTHLY	BI-MONTHLY	SEMI-ANNUAL	ANNUAL	AS NEEDED
<b>Vacuum Pump</b>								
a.)Check & Clean Filter								
b.)Check For Water Leaks								
c.)Run Cleaner through Lines								
<b>X-Ray Units, Intra-Oral</b>								
a.)Check Arms for Drift								
b.)Check Tubehead for Oil Leaks								
c.)Verify Audible Alert								
d.)Calibrate Every Two Years								
e.)Check Lead Aprons								
<b>X-Ray Unit, Panoramic</b>								
a.)Check Audible Alert								
b.)Check Cassettes								
c.)Check Cassette Holder								
d.)Calibrate Every Two Years								
e.)Check Lead Aprons								
<b>Whisper Jet, Air Abrasion</b>								
a.)								