

EQUIPMENT MAINTENANCE

DESCRIPTION OF EQUIPMENT								
Automatic Film Processor								
a.) Check Chemistry	■							
b.) Change Chemistry				■	■			
c.) Clean Racks d.) Clean Entire Unit		■	■			■		
Chairs								
a.) Clean		■	■					
Compressor								
a.) Check Moisture Indicator			■	■				
b.) Check Air Intake Filters				■	■			
c.) Regenerate Dessicant							■	■
d.) Check Oil Level		■	■					
e.) Change Oil							■	■
Dental Operating Light								
a.) Clean Reflector							■	■
b.) Check and Clean Sheild		■	■					
c.) Check for Drift				■	■			
Nitrous System								
a.) Check Pressure Overnight						■	■	
b.) Turn System On/Off	■	■						
c.) Check Rubber Goods for Wear					■	■		
Sterilizer, Pelton								
a.) Activate Safety Valve b.) Clean Chamber c.) Spore Test d.) Check Door Seals					■	■	■	
Sterilizer, Statim								
a.) Lubricate Cassette seal			■	■				
b.) Apply Stat Dry to Cassette								
c.) Replace Air Filter						■	■	
d.) Empty Used Water Container							■	■
Dental Unit								
a.) Clean Out Vacuum Trap	■	■						
b.) Clean Vacuum Lines			■	■				
c.) Flush Handpiece Tubings		■	■					

